

HYDRANAL SMART REAGENTS FOR METROHM OMNIS KF TITRATORS



Hydranal™ reagents and Metrohm instruments are industry leaders & collaborators in Karl Fischer titration. We are now excited to exclusively launch a special Hydranal product line for Metrohm's OMNIS KF titrators.

Our Hydranal Smart Reagents bottles are embedded with an RFID chip sticker that captures essential information for every titration, including article code, product name, lot number, shelf life and numeric concentration. Connect Hydranal reagents to your Metrohm OMNIS KF titrator in a faster, safer and better manner. Simply peel off the add-on RFID sticker from the Hydranal bottle and attach to the 6.01601.000 Bottle Cap multi-use in 3 easy steps!



FASTER

Automated data transfer

SAFER

Reduced chemical contact

BETTER


End-to-end traceability

Honeywell


 **Metrohm**

INSTRUCTIONS FOR RFID STICKER USE


RFID sticker to be used on METROHM OMNIS Bottle cap multi-use



1 Pull out RFID sticker and peel off.
Note: Do NOT bend the sticker.



2 Apply RFID sticker on OMNIS Bottle cap multi-use
Note: Remove any old stickers from cap before applying.



3 Secure OMNIS Bottle cap multi-use on Hydranal bottle and connect with the OMNIS Liquid Adapter.
Note: For accurate results always perform an initial titer determination every time you connect a new reagent

HYDRANAL SMART REAGENTS FOR METROHM OMNIS KF TITRATORS

PRODUCT NAME	PRODUCT NUMBER	DESCRIPTION	PACKAGING
HYDRANAL NEXTGEN Composite 5 OM	34805-1L-OM	One-component reagent, titer ~5 mg/mL	1 L
HYDRANAL NEXTGEN Titrant 5 OM	34801-1L-OM	Two-component reagent, titer ~5 mg/mL	1 L

NOTE:

- Hydranal Smart Reagents are compatible with Metrohm OMNIS titrators with software version 2.10 or higher. You can easily download the software update if required.
- There is no change in the chemical composition of Hydranal reagents when used as Hydranal Smart Reagents. Customers wishing to upgrade to Hydranal Smart Reagents can do so without any corresponding change in their method. An official confirmation can be found on our website.



PT. INDOFA UTAMA MULTI CORE

IMPORTER – AGENT - DISTRIBUTOR

Laboratory & Scientific Equipments • Chemicals • Glassware

SURABAYA HEAD OFFICE: Jl. Mawar No.27 - 29 Surabaya – 60262 Telp: (031) 5329080

Fax: (031) 5329077 ; (031) 5329078 e-mail: iumcsby@indofa.com

BRANCH OFFICE: PT. INDOFA UTAMA MULTI CORE
Jl. Kemanggisan Raya 125,
Palmerah Jakarta Barat- 11480 Telp: (021) 53666156
Fax: (021) 5484120 e-mail: jakarta@indofa.com

REPRESENTATIVE OFFICE: CV. JURUS MAJU
Jl. Seteran Serut No.4 Semarang – 50134
Tel. (024)3555693
Fax. (024) 3552543 email: sc@jursmaju.com

All statements and information provided herein are believed to be accurate and reliable, but are presented without guarantee, warranty or responsibility of any kind, express or implied. Statements or suggestions concerning possible use of our products are made without representation or warranty that any such use is free of patent infringement, and are not recommendations to infringe any patent. The user should not assume that all safety measures are indicated herein, or that other measures may not be required. User assumes all liability for use of the information and results obtained.



RESPONSIBLE CARE[®]
OUR COMMITMENT TO SUSTAINABILITY

For more information, visit

hydranal-honeywell.com
metrohm.com/omnis-kf

Hydranal Center of Excellence

Tel. +49 5137 999 353
Fax: +49 5137 999 698
hydranal@honeywell.com

6116 | AM | 06/20
© 2020 Honeywell International Inc.

Honeywell

Metrohm