

JP Recirculating Water Pump for Fumes Aspiration



GENERAL FEATURES

Construction material:	ABS
Air flow max:	35 l/min
Flow regulation:	from 0 to 35 l/min
Residual pressure with water temperature of 15°C:	35 mmHg
Power:	160 W
Power supply:	115 V or 230 V / 50-60 Hz
Weight:	8,4 Kg (18.5 lb)
Dimension (WxHxD):	250x400x370 mm (9.8x15.7x14.6 in)

BENEFITS

Automatism and optimization of the sucking power

Strongly reduced water consumption in comparison to water jet pumps connected to tap

Possibility of avoiding emission into the atmosphere of gaseous pollutants or polluted waters into drainages

Independent from tap water contingencies like low water availability or low pressure