

SMS Scrubber



GENERAL FEATURES

Construction material:	Epoxy painted stainless steel structure
Capacious containers for the neutralization solution and for the condensate collection	
Weight:	3,5 Kg (7.7 lb)
Dimension (WxHxD):	300x500x190 mm (11.8x19.7x7.5 in)

BENEFITS

- A highly efficient neutralization, due to the optimized gas-liquid contact
- Used for both basic and acid neutralization
- The apparatus is safe and practical to use
- Simple disposal of condensates
- Simple control and replacement of neutralizing solution