



SZ AIR TECH



Laminar Air Flow Vertical and Horizontal Flow



Desk-Top Model



VD-650 Vertical	Model	HD-650 Horizontal
Class 100(FED 209E)	Flow Direction	Class 100(FED 209E)
Adjustable 0.3~0.5m/s	Cleanness	Adjustable 0.3~0.5m/s
≤62dB(A)	Air Velocity	≤62dB(A)
≤5μm	Noise	≤5μm
≥300Lux	Vibration	≥300Lux
220V, 50Hz, 1Phase	Illumination	220V, 50Hz, 1Phase
300Watts	Power Source	300Watts
60	Power Consumption	55
615x495x500mm	Weight (kg)	615x400x535mm
650x535x835mm	Working Zone (wxdxh)	650x680x625mm
1x8W	Overall Dimensions	1x8W
1pc – 610x450x50mm	Number of Light	1pc – 610x490x50mm
Single	Hepa Filter	Single
	Number of Operator	



Vertical Flow Model



SW-CJ-2F/1F



SW-CJ-2BU/1BU

Horizontal Flow Model

SW-CJ-1F	SW-CJ-2F	Model	SW-CJ-1BU	SW-CJ-1CU
Class 100(FED 209E)	Class 100(FED 209E)	Cleanness	Class 100(FED 209E)	Class 100(FED 209E)
≤0.5/Ø90mm	≤0.5/Ø90mm	Number of Colony	≤0.5/Ø90mm	≤0.5/Ø90mm
≥0.3m/s	≥0.3m/s	Air Velocity	≥0.3m/s	≥0.3m/s
≤62dB(A)	≤62dB(A)	Noise	≤62dB(A)	≤62dB(A)
≤3μm	≤5μm	Vibration	≤3μm	≤4μm
≥300Lux	≥300Lux	Illumination	≥300Lux	≥300Lux
220V, 50Hz, 1Phase	220V, 50Hz, 1Phase	Power Source	220V, 50Hz, 1Phase	220V, 50Hz, 1Phase
1100Watts	1500Watts	Power Consumption	400Watts	650Watts
150	300	Weight (kg)	110	200
870x690x520mm	1360x690x520mm	Working Zone (wxdxh)	820x480x600mm	1680x480x600mm
1000x750x1600mm	1490x750x1600mm	Overall Dimensions	900x700x1450mm	1760x700x1450mm
1x15W & 1x15W	2x15W & 2x15W	Number of Light & UV	1x20W & 1x20W	1x40W & 1x40W
1pc – 820x600x50mm	2pc – 610x610x50mm	Hepa Filter	1pc – 820x600x50mm	2pc – 820x600x50mm
Single	Double	Number of Operator	Single	Double



VS-series

Detachable Model



VT-series



HS-series

Detachable Model



HT-series

VS-840K-U	VS-1300L-U	VT-840-U	VT-1300-U	Model	HS-840-U	HS-1300-U	HT-840-U	HT-1300-U
ISO Class 5 Class 100(FED 209E)				Clean-ness	ISO Class 5 Class 100(FED 209E)			
≤0.5/Ø90mm				No of Colony	≤0.5/Ø90mm			
≥0.3m/s				Air Velocity	≥0.3m/s			
≤62dB(A)				Noise Level	≤62dB(A)			
≤4μm				Vibra-tion	≤3μm			
≥300Lux				Illumi-nation	≥300Lux			
220V, 50Hz, 1Phase				Power Source	220V, 50Hz, 1Phase			
400Watts	650Watts	400Watts	650Watts	Power Con-sum-ption	400Watts	650Watts	400Watts	650Watts
230	420	230	420	Weight (kg)	160	350	160	350
660x650x570 mm	1120x650x720 mm	660x650x570 mm	1120x650x720 mm	Working Dim (wxdxh)	720x650x570 mm	1180x650x570 mm	720x650x570 mm	1180x650x570 mm
840x825x1610 mm	1300x825x1760 mm	840x825x1610 mm	1300x825x1760 mm	Overall Dimensions	840x825x1440 mm	1300x825x1440 mm	840x825x1440 mm	1300x825x1440 mm
1x20W & 1x20W	1x30W & 1x30W	1x20W & 1x20W	1x30W & 1x30W	No of Light & UV	1x20W & 1x20W	1x30W & 1x30W	1x20W & 1x20W	1x30W & 1x30W
1x 760x610x50mm	2x 610x610x50mm	1x 760x610x50mm	2x 610x610x50mm	Hepa Filter	1x 760x610x50mm	2x 610x610x50mm	1x 760x610x50mm	2x 610x610x50 mm
Single	Double	Single	Double	No of Opera-tor	Single	Double	Single	Double

Untuk informasi lebih lanjut, harap hubungi:

